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The product is intended to be used in human beings for the purpose of treating or alleviating a disease, as well as for the treatment, alleviation, or compensation of an injury or disability, in accordance with the current Spanish regulations on medical devices.

According to the UNE-EN 60601-2-52 standard for medical beds, it may be intended **for two different application environments (AE)**:

- EA3.- long-term care provided in a medical area where medical supervision is required and monitoring is provided if necessary, and where the bed is used in medical procedures that may be provided to help maintain or improve the patient's condition.
- EA4.- care provided in a domestic area where the bed is used to alleviate or compensate for an injury, disability, or disease.

The product is reusable. The expected service life is estimated at 10 years under normal conditions of use. Actual service life may vary depending on usage intensity and compliance with the intended use. The stated service life applies provided the product is used on a daily basis and in accordance with the safety, cleaning, storage, and maintenance requirements specified in the Instructions for Use (IFU)

The bed is considered a Class I medical device in accordance with European Regulation (EU) 2017/745 on Medical Devices (MDR).

Reference Standards:

*** UNE-EN 60601-2-52** "Medical electrical equipment. Particular requirements for basic safety and essential performance of medical beds"

**For continuous product improvement, Technimoem Care, S.L. reserves the right to modify technical specifications and/or product images without prior notice.*

TECHNICAL SPECIFICATIONS		
EXTERNAL DIMENSIONS AND WEIGHTS		For 90 × 190 cm mattress
Bed weight:		37 Kg
External bed dimensions:		90x196 cm
Bed with 4-bar side rail set:		102x196 cm
Bed with 4-bar side rail set and Corera headboard and footboard:		102x202 cm
OTHER MAXIMUM PERMISSIBLE LOADS (ESTIMATED)	150 kg. Maximum user weight	
	25 kg. Mattress + bedding weight	
	15 kg. Side rails weight (approx.)	
	35 kg. Weight of remaining accessories (lifting pole, headboard, footboard, etc.)	
LIFTING FRAME TYPES AND WEIGHTS (versión plus)		
Weight of lifting frame with independent brake:		36,00 kg
Weight of lifting frame with centralized brake:		40,40 kg
FRAME-TO-FLOOR HEIGHT ADJUSTMENT RANGE		From 38 a 80 cm
Wheels Ø125 mm		
LEG TYPES (Ø50 mm tube)	Mattress platform height	Weight (per set)
Adjustable legs from 35 to 46 cm	From 41 to 52 cm	3,00 Kg
Fixed braced legs 35cm	41 cm	4,00 Kg
Adjustable braced legs from 35 to 50cm	From 41 to 56 cm	6,20 Kg
Fixed braced legs with Ø125 wheels, 35 cm	41 cm	7,00 Kg
Adjustable braced legs with wheels from 35 to 45 cm	From 41 to 51 cm	8,20 Kg

BED BASE

External frame made of 50x20x1.5 mm steel tube and internal joints made of 25x25x1.5 mm steel tube, curved at the corners (manufactured according to UNE 10305-3 standard).

30 mm safety clearance around the entire perimeter of the bed between the joints and the external frame to prevent entrapment.

Epoxy-polyester paint coating.

Equipped with square fittings in the headboard area for attaching an IV pole holder and/or lifting pole.

Option to install side safety rails (consult models).

Option to install wooden headboard and/or footboard (consult models).

Metal mattress retainer arch at the foot of the internal joint and injected thermoplastic mattress retainers (PP with fiberglass) on the sides.

- Wood slat bed base option: 17 curved steamed beech slats (53x8 mm), coated with melamine paper (3 of them double with tension adjusters for lumbar reinforcement). Slat holders and adjusters made of rigid polyethylene in light grey.
- Injected thermoplastic slat bed base option: (polypropylene with fiberglass PPFV).


4-SECTION ELECTRIC ADJUSTABLE BED (3 ARTICULATIONS):

Angle A: Backrest adjustable from 0° to 70° relative to the horizontal.

Angle B: Thigh section adjustable from 0° to 20° relative to the horizontal.

Angle C: Calf section adjustable from 0° to 10° relative to the horizontal.

The thigh and calf sections are adjusted simultaneously.

ELECTRICAL EQUIPMENT

Both models come standard with the same lifting carriage motor as the Plus version, allowing for easy transformation of the bed from a leg-based bed to a trundle bed.

All functions are controlled with a single remote.

Bed articulation: Via a minimum 4500+4500N electric motor with a wired remote. 230V input. 24V output (with an internal 50W transformer). 581mm distance between axles. 87/69mm stroke. IPx4 protection rating. Class II.

Supplied with one or two 9V batteries in case it is necessary to return to the horizontal or resting position in the event of a power outage.

For NULES PLUS

Trunk movement: Via a linear actuator with a minimum of 6000N, with an installation distance of 500mm and a stroke of 250mm, connected to the electric motor.

UNDER-BED LIGHT: Easily illuminates the bed's perimeter using the remote.

REST POSITION: Place the bed in the rest position by touching a single button. The bed lies flat and, in the case of the PLUS version, lowers in height simultaneously.

TYPES OF LEGS

- Adjustable – Plastic foot

Ø50 mm tube with leveler.

Height: 35 to 46cm

Also available in 30 to 41 cm and 25 to 36 cm options.

- Fixed braced legs

Ø50 mm tube joined with 30x15x1.5 mm tube.

Height: 35 cm.

- Braced adjustable legs with screw

Ø50 mm tube joined with 30x15x1.5 mm tube.

Height: 35 to 50 cm.


- Fixed braced legs with castors

Ø50 mm tube joined with 30x15x1.5 mm tube.

Castors with brake Ø125 mm.

Fixed height: 35 cm.

- Braced adjustable legs with castors

Ø50 mm tube joined with 30x15x1.5 mm tube.

Castors with brake Ø125 mm.

Height: 35 to 45 cm.

ELEVATION SYSTEM, for the NULES PLUS model

External frame made of 40x40x2 mm and 40x30x1.5 mm steel tubing.

Internal joints made of 50x25x2 mm steel tubing (cross-type design).

Bed connection arms made of 40x30x1.5 mm steel tubing.

Side crossbars made of 50x20x1.5 mm, 25x25x1.5 mm, and 40x30x1.5 mm tubing.

All tubes are manufactured according to UNE 10305-3 standard.

Epoxy-polyester paint coating.

Clearance under the bed: 160 mm, allowing the use of a patient lifting hoist.

4 independent braked castors, Ø125 mm.

Optional frame with centralized braking system operated by a stainless steel lever and one directional castor.

Height adjustment range: from 38 cm to 80 cm.

LIFTING POLE (optional)

Structure made of Ø33.7x3.2 mm steel tubing coated with epoxy-polyester paint. Includes a plastic handle and polyamide fabric straps. Maximum safe working load: 75 kg.

Weight: 4,4 Kg


IV POLE (optional)

Structure made of Ø33.7x3.2 mm steel tubing coated with epoxy-polyester paint.

Equipped with two hooks for hanging two bags.

Weight: 3 kg.

CLEANING AND MAINTENANCE

The bed has been designed to allow easy cleaning and optimal disinfection. Follow these instructions to prevent deterioration that could compromise the proper functioning of the bed and void the product warranty.

- Ensure that the bed is immobilized and disconnected from the mains power supply.
- Keep castors and brakes free from dirt and soap residues.
- Never wash the bed with large amounts of water, high-pressure jets, or in an automatic washing tunnel.
- Do not use water at temperatures exceeding 60°C.
- Avoid excessive exposure of connection plugs to water.
- Carefully dry the bed before returning it to service.
- Use clean single-use or recyclable cloths.
- Use a detergent solution or a disinfectant solution containing 70% diluted alcohol.
- Use non-abrasive products (pH between 5 and 9).

BEFORE ANY USE OF THE BED In addition to the cleaning measures described above, any bed components showing defects must be replaced prior to use. Contact the distributor to obtain the necessary spare parts.

Maintenance operations may only be carried out by qualified and authorized technical personnel. Before performing any repair, maintenance, cleaning, or disinfection procedures:

- Disable the electrical functions.
- Disconnect the bed from the main power supply if no electrical operation is required.
- Immobilize the bed and take all necessary precautions to prevent movement.
- Never open or puncture an electric motor.
- Contact the supplier for any specific maintenance issues.
- Refer to the manufacturer or visit www.tecnimoem.com for component replacement instructions.

An annual periodic inspection of the bed and its accessories is recommended to ensure proper operating condition, extend the product's service life, and contribute to sustainability objectives. Particular attention should be paid to:

- Bed lifting and articulation system
- Mains power connection cable
- Hand control unit
- Overall condition of the bed
- Headboard and footboard
- Side rails
- Castors
- Brake pedals
- Emergency batteries
- Underbed light

Carry out all necessary repairs and replace components where required.

SPARE PARTS LIST

The periodicity and frequency of technical inspections must be adapted to the condition of the bed and its specific conditions of use, as certain components may require periodic replacement.

Upon request, all plastic components, electrical equipment, slats, castors, and accessories may be replaced.

To ensure the declared service life of the bed, TECNIMOEM recommends replacing the following components:

Description	Frequency
Wheels	Every 10 years
Bed motor	Every 5-7 years
Lift motor	Every 5-7 years
Hand control	Every 5-7 years
Under-bed light	Every 5-7 years
Plastic slat	Every 3-5 years
Plastic support block	Every 3-5 years
Mattress retainer	Every 3-5 years
Headboard/footboard hardware	Every 5 years
Wooden headboard/footboard	Every 10 years
Metal frame	Every 10 years

TRANSPORT AND STORAGE

All necessary precautions must be taken to ensure the safe transport of the bed and its accessories, preventing impacts and exposure to dust.

During transport and/or prolonged storage, the bed must be kept in the lowest position, with all functions deactivated, castors braked, secured with straps, protected against water ingress and humidity (< 75%), and stored at a temperature between -20°C and +50°C.

If the bed has been stored for an extended period (more than 4 months), its operation must be inspected by qualified personnel prior to being returned to service.

ENVIRONMENT AND END OF LIFESPAN

At the end of the product's lifespan, it must be taken to a recycling centre. This helps contribute to environmental care.

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Tecnimoem Care develops, manufactures, and distributes beds, bed bases, hoists, chairs, and accessories for geriatric and domestic use.

In a world where the quality of life reaching old age is the most valued asset, **Tecnimoem Care** puts all its effort into continuous research and development of elements for rest that people who due to physical problems, age or simply who value ergonomics and adaptability, choose our products.

Tecnimoem Care, in its innovation effort, incorporates the best materials on the market with perfect finishes that exceed all quality standards and are what differentiate us and make us be chosen by our customers.

SRN: ES-MF-000019199
Según (EU) 2017/745 MDR

TECNIMOEM CARE S.L.
CTRA. Logroño-Mendavia Km.5
Tfno. +34 948 64 62 13
31230 VIANA (Navarra) – ESPAÑA (SPAIN)
tecnimoem@tecnimoem.com
www.tecnimoem.com